

**Particle Size Analysis  
in turbid suspensions and emulsions  
from 1 nm to 10.000 nm**

**NANOPHOX**

**HST: Cluster Abel 1689 Galaxy**




# Introduction

★ Particles in the **nanometre range** offer numerous new properties:

- ✓ Scratchproof surfaces and coatings
- ✓ Conductive glass
- ✓ Non corrosive steel
- ✓ Optical storage capacity of plastics
- ✓ Bullet proof synthetic fibres

★ Can be designed e.g. in composite products, which contain large amounts of transparent fumed silica in the nanometre range as filler in a plastic matrix




Example Aerosil Typ „0.2 – 2  $\mu\text{m}$ “  
about. 50 Weight. %




# Size Scales

★ The new properties of these fillers are related to their particle size




# Material Characteristics

- ★ Material characteristics depend on
  - ✓ The **particle size distributions**
  - ✓ The **stability** of such distributions
- ★ Nano particles having extremely large surfaces tend to




# Photon Correlation Spectroscopy (PCS)

## Interaction of laser light and particles



# Interaction of laser light and particles in the nanometre range

## ★ Brownian molecular motion



- ★ The changing short-range order of the particles yields a fluctuating intensity distribution
- ★ The intensity related to the scattering angle  $\Theta$  oscillates with time  $t$



★ Brownian molecular motion

- ✓ The phenomenon was correctly assumed by W. Ramsay in 1876 and proven by A. Einstein and M. von Smoluchowski in 1905/06
- ✓ „Stokes Einstein“ relation applies:

$$D(x) = k_B T / 3\pi\eta x$$

D	diffusion constant
$k_B$	Boltzmann-constant
T	absolute temperature
$\eta$	dynamic viscosity of surrounding medium
x	particle diameter

applying the auto-correlation function  $G(\tau)$

$$G(\tau) = \langle I^s(\theta, 0) I^s(\theta, \tau) \rangle = \langle I \rangle^2 (1 + \exp(-2 D(x) q^2 \tau))$$

$$\langle I \rangle^2 (1 + \exp(-2 q^2 k_B T \tau / 3\pi\eta x))$$

↪ decay constant  $\propto$  diffusion  $D(x) = \frac{k_B T}{3\pi\eta x}$



↪ decay constant  $\propto$  1/particle size x



★ Determination of particle size

- ★ Calculation of auto correlation function  $G(\tau)$
- ★ Evaluation of particle size from slope  $\ln(G(\tau)-1)$

intensity  $I(t, \Theta)$





## Limitations of PCS

- ★ Valid only for **non-interacting, spherical** particles
- ★ **Singularly scattered** laser light
- ★ **Multiple scattering** in turbid suspensions of high concentrations **distorts** the result
  - ✗ Depends on **concentration**
  - ✗ And cuvette position
- ★ To avoid multiple scattering PCS has to be operated at **extremely low concentrations**:
  - ✗ Low scattering intensities
    - ↪ **Poor** signal-/noise ratio
    - ↪ **Long** analysis times
  - ✗ Dilution not only delicate → often even **non-permissible**
  - ✗ Contaminations → **dust bastardises** the result
  - ✗ **Stability** of sample ?




Improvement of principle: Elimination of influence of **multiple scattering**



# Photon Cross Correlation Spectroscopy (PCCS)

## Principle



✓ Two laser beams cross over in the sample container and generate two similar signal patterns (speckles)


✓ Identical parts can be selected by correlation of the intensities of two equivalent signal spots



# Result

- ★ Key principle: 3D-cross correlation
- ✓ Selection of singularly scattered light

2 measured variables:



1. Amplitude  $\propto$  singularly scattered part

✓ Mutation of sample  $\rightarrow$  change of amplitude

↪ Potential for stability analysis


2. Slope  $\propto$  1/particle size  $\times$

✓ PCCS provides exact information for turbid samples




# Elimination of multiple scattering portion:


★ Latex in H<sub>2</sub>O, 107 ± 10 nm (TEM)



Transparent: Transmission 99.7 %  
 Turbid: Transmission 0.7 %  
 Both: (1x1 cm<sup>2</sup>-cuvette)




Copyright©2005, IFAM Bremen, Dr.M.Kleemeier




Copyright©2000, ALV-GmbH, Langen Application Report #004



# Independence on Cuvette Position and Concentration:




Latex 120 nm,  
0.2 weight %, 3 different cuvette positions



Latex 120 nm,  
4 concentrations: 2; 0.4; 0.2; 0.1 weight %



# Concurrent Analysis of Stability and Particle Size




Sedimentation of a turbid quartz suspension over many hours




# Evaluation

## Standard evaluation (2<sup>nd</sup> Cumulant)

- ★ Presumes and requires mono modularity and yields:
  - ✓ Median **diameter** of the distribution
  - ✓ Indication as to the **width** of distribution (width  $\pm$  %)
  - ✗ No information regarding the **internal structure** of the total distribution




mean diameter 195 nm, width 10,2 %



mean diameter: 195 nm  
width:  $\pm$  10.2 %




# NNLS (Non Negative Least Square):



✓ Selection of the best fitting window determines the quality of the evaluation

✓ Indication of structure of the distribution possible

✓ For selection of the properly detected area (best fit) reliable information for the single modes





# WINDOX 5 Operation Control and Evaluation Software

**NANOPHOX**

File View Measurement Database Application Extras Tools Help

Measurement: 2004-11-05 14:41:22.2140 0102 P

Product: Latex  
Liquid: Water  
Measuring cond.: 60s 150kcps 20x  
Evaluation mode: NNLS 10..1000nm:128(log) 0.0010..10.0000ms 2.0%

9

Output

**Start of evaluation:** Evaluate current measurement  
 Immediately after measurement

**Report:**  
Screen:  Measurement  Statistics  
Printer:  Measurement  Statistics

**Distribution-diagram:**  Measurement  Statistics

**Q(t) diagram:**  Measurement

**Correlation-function:**  Measurement

WINDOX 5.1.2.2

**Nanophox signal test**

Mode View

Static light scattering:  
Single scatt. ratio 0 %

**Count rates:**  
∅ A (total) 122.0 kcps  
∅ B (total) 166.5 kcps  
∅ A+B (moving) 144.3 kcps

**Device status:**  
Probe temperature 25.00 °C  
Desired temperature 25.00 °C  
Laser intensity 19 %  
Motor position X 12.0 mm  
Motor position Y 8.0 mm

**Particle size:**  
Particle diameter 124.53 nm  
Width of distribution 0.00 %

**Decay constant:**  
Decay constant 0.00 nm

A B A+B Ex RMS

TRACE

Count rate [kcps]

Count rate channel A Count rate channel B

t[s]

Cor

Correlation function


Correlation

t[ms]



# NANOPHOX™

PCCS – technology made by Sympatec GmbH



## Technical Data NANOPHOX™

Measuring range:	1 nm – 10.000 nm
Principle:	Dynamic light scattering
Scattering angle:	90 degree
Evaluation:	3-D cross correlation
Analysis volume:	ca. 0.3 ml
Sample container:	Cuvette 10x10 mm <sup>2</sup>
Cuvette position:	Automatic or manual
Light source:	HeNe-laser, 10 mW max.
Wave length:	632.8 nm
Class of laser:	3B
Class of instrument:	1
Optics:	No adjustment requested
Temperature range:	15 – 40°C , typ. 22 – 25°C recommended, via software
Temperature stability:	0.05°C ( $\cong$ 0.1 % precision during evaluation of particle size)
Precision of temperature control:	0.1°C $\cong$ 0.2 % precision of particle size



# NANOPHOX Conclusions

★ Photon Cross Correlation Spectroscopy (PCCS) provides:

- ✓ Concurrent analysis of **particle size** from **1 nm** to **10 μm** and **stability** of sample
- ✓ For T=const. only **η** (dynamic viscosity of liquid) requested (refractive index **n** for mass and volume proportions only)
- ✓ Cross correlation completely **eliminates** the **multiple scattering allotment**
  - ↪ Analysis with **turbid** and **coloured** samples
  - ↪ Independence of results on **concentration** and **cuvette position**
  - ↪ No adaptation of evaluation to multiple scattering necessary
- ✓ Analysis for **high concentrations**
  - ↪ Normally **no dilution** required
  - ↪ No special request for cleanliness (as for PCS)
  - ↪ No special glass cuvette necessary)
  - ↪ **Rapid, statistically secure analysis** due to high counting rate
  - ↪ Easiest **handling**

