

**Particle Size Analysis and Shape Evaluation
with high speed image analysis
and efficient dispersion**


QICPIC

HST: HELIX Nebula



Introduction

- ★ For characterisation of *real* particles *particle shape* in addition to the *particle size* gains more and more importance
 - ☆ Image Analysis (IA) can give valuable services
- ★ *Static image analysis* (non-moving particles, e.g. on a microscope slide)
 - ✓ Depth ϵ of the sharpness is defined
 - ✓ High resolution for small particles
 - ✓ Well established and standardised in ISO 13322-1 (Quimet, Kodak)
 - ☆ Small amount of data
 - ✗ Particles are *orientated* by the base
 - ✗ Problem with *overlapping particles*
 - ↪ High effort for software correction
 - ✗ *Tiny sample size* → sampling problem
 - ↪ Very poor statistics



Dynamic Image Analysis (Particles moving)

★ A flow of moving particles is captured


- ✓ *Enlarged sample size*
- ✓ *Arbitrary orientation*
- ✓ *Reduction of overlapping particles*

★ Several companies provide systems

- ☆ Reflection or transmission technology
- ☆ Wet dispersion or free fall
- ☆ Matrix or line scan cameras

★ Common to all systems is

- ✗ *Dry products: Free fall* → for well flowing bulk materials only, no dispersion
- ✗ *Wet dispersion* → *smallest sample sizes* and *slow particles*
- ✗ *Low frame rates* → small number of particles → *bad statistics*
- ✓ All system types will be integrated in the coming standard ISO 13322-2



Optimum



Dynamic image analysis with *powerful dispersion* and *good*




New Approach (1): Combination with Dispersers

- ★ Effective **dispersion** must be the **prerequisite** of image analysis
- ★ This adds **energy** to the particles → creates fast particles
→ which have to be **imaged clearly** within **exposure time**

Dry disperser RODOS™

- ★ Combines particle-to-particle, particle-to-wall and centrifugal forces caused by velocity gradients for effective dispersion down to **$0.1 \mu\text{m}$**
 - ★ Particle velocity in aerosol beam: Up to **100 m/s**
 - ★ Flash lamps with $\tau = 100 \mu\text{s}$: Motion blur up to **10 mm (!)**
 - ★ Best flash lamps with $\tau = 1 \mu\text{s}$: Motion blur of up to **$100 \mu\text{m}$**
- ↳ Required exposure time:
 $\tau < 10 \text{ ns}$ for motion blur of **$< 1 \mu\text{m}$**



New Approach (2): High Frame Rates

- ★ In order to obtain a *representative sample* and a *good statistical reliability* a **large number** of particles should be acquired in **short time**

Example 1 1 % precision requires 10.000 particles within one size class
 Within 10 particles/image → 1000 images or 40 s measuring time
 at 25 images/s per class. → for 30 equally populated classes:
 The measuring time is about *1200 s = 20 minutes*

↪ *Not acceptable*

Example 2 As above, but 500 images/s, and because of good dispersion
 20 particles/image → for 30 equally populated classes:
 The measuring time is only *30 s* !


↪ *OK*

↪ *High frame rate required* (>> 25 images/s)
 ↪ *Extreme data volumes* (500 Mbyte/s at 500 images/s, 1 Mpixel, 8 bit)




New Approach (3): Special Optical Set-up

- ★ Even a *perfect parallel beam* will **not** produce a clear shadow of a small particle
 - ↪ The intensity pattern depends on the distance between sensor and particle



- ★ Without imaging optics
 - ✗ Sensor must be *very close* to the particles
 - ✗ Image quality *depends on particle size* → high calculation effort
 - ✗ Sensor must be scaled to the particle size


Standard Imaging Lens



- ★ A lens transfers light from one point at the *object plane* to one point in the corresponding *image plane*
- ★ At image plane 2' the image of arrow 1 is blurred
 - ↳ Only objects located at a single object plane are in focus on a flat sensor
- ★ The size of *blurred spot* depends on *distance of the object to the object plane* and the *angle of rays* θ which passes the lens
- ✗ **Magnification** depends on **object position** (perspective image)




Special Imaging Lens





- ★ Only ray fans *parallel to the optical axis* are used for the image
- ★ The aperture stop controls the angle of rays θ (aperture angle)
- ✓ Image size does not depend on object position
 - ↪ Image size *less sensitive to defocus*
 - ↪ Even *transparent particles* show *high contrast*, as deflected light is not imaged
 - ↪ Calculation effort is *remarkably reduced*

transparent glass beads



Illumination




- ★ Front illumination
 - ✓ Very *realistic images incl. surface structure*
 - ✗ Reflection depends on *particle material and shape*
 - ✗ Difficult and error-prone image analysis

- ★ Illumination box
 - ✗ Light can still be *reflected by particle borders*
 - ✗ *Small aperture stop* is required to reduce reflection
 - ↪ Very much light is required
 - ↪ *Limited imaging* of small particles




Illumination in Transmission



- ✓ The *smallest aperture angles* can be used → **highest contrast**
- ✓ The aperture stop blocks only *stray light* and large angle *diffracted light* from the particles
- ✓ The aperture stop has *no influence* on the brightness of the background
→ only **small light intensities** are required
- ✓ The *depths of field* is limited by diffraction only → **wide working distance**





★ Pulsed light source

- ✓ Pulse duration: **< 1 ns**
- ✓ Output power: $\approx 0,15$ nJ/pulse
- ✓ Light colour: visible (≈ 532 nm)
- ✓ Repetition rate: **0 to 500 Hz**, adjustable

★ Beam expansion unit

- ✓ 35 mm Ø, 16 mm Ø, 7 mm Ø


★ Imaging objectives

- (up to 5 on a carousel)
- ✓ 2:1, 1:1, 1:3, (1:10) ...

★ Camera

- ✓ CMOS, 1024 x 1204 pixel,
- ✓ Square pixel, $10 \mu\text{m}^2$
- ✓ Up to **500 fps**






QICPIC & RODOS/L




QICPIC & RODOS/L



Measuring zone of
dry disperser RODOS/L





QICPIC & GRADIS/L


Gravity disperser



Side view



QICPIC & OASIS/L





dry & wet
disperser





QICPIC & LIXELL

wet disperser



QICPIC Front View

With open camera



✓ Simplifies cleaning and exchange of dispersion unit



QICPIC Sensor Control

QICPIC

File View Measurement Database Application Extras Tools Help

Measurement: 2003-11-25 16:40:56.7500 Q0103

Product: Silicon carbide P50

Measuring range: 1.00 63.0 mm

Trigger condition: Top 0.2%

Disperser: GRADIS GRADIS&VIBRI_25%

Sample: Sample P50'01 - 03

Trigger Condition

Normal measurement

Start 0.000s after opt. concentr. \geq 0.2%

Valid

Stop after 10.000s opt. concentr. \leq 0.2%

or after 20.000s time. Trigger timeout: 60s

WINDOX 5.1.1.2

WINDOX5 Database browser <QICPIC Database 2>

All Products All Ranges All Dispersers All Dates

Date	Time	Product name	Measuring range	Dispersing system	Evaluation	Parameter1	Review
2003-11-25	13:53:00.9680	ball bearings	1.00 63.0 mm	GRADIS	USER_1	calibration 01	
2003-11-25	14:33:26.3590	ball bearings	1.00 63.0 mm	GRADIS	USER_1	calibration 02	
2003-11-25	14:47:52.7960	ball bearings	1.00 63.0 mm	GRADIS	USER_1	calibration 03	
2003-11-25	14:50:18.3590	ball bearings	1.00 63.0 mm	GRADIS	USER_1	calibration 04	
2003-11-25	14:56:50.2180	ball bearings	2.00 63.0 mm	GRADIS	USER_1	calibration 05	
2003-11-25	16:33:45.9840	Silicon carbide P...	1.00 63.0 mm	GRADIS	EQPC	Sample P50'01 ...	
2003-11-25	16:36:34.8750	Silicon carbide P...	1.00 63.0 mm	GRADIS	EQPC	Sample P50'01 ...	
2003-11-25	16:40:56.7500	Silicon carbide P...	1.00 63.0 mm	GRADIS	EQPC	Sample P50'01 ...	

Total:18 Selected:1 Current:8


★ Every measurement is characterised with:

- ☆ Unique *primary key*
- ☆ Product
- ☆ Measuring range
- ☆ Trigger condition
- ☆ Disperser
- ☆ User specific parameters




Signal Test

- ✓ *On-line monitoring* during measurements (10 fps)
Video player function in database mode
 - ☆ Scaled and non-scaled
- ✓ Selection of a particle directly:
Displays size and shape parameters
- ✓ *Export* via clipboard or video streaming file (AVI)



Calculation Modes




- ✓ Circle of equal projection area
- ✓ Feret Diameter
 - * maximal
 - * minimal
 - * averaged
 - * 90° to maximum
 - * 90° to minimum
- ✓ Minimal area of enclosing rectangle
 - * length
 - * width
- ✓ Chord length
 - * vertical
 - * horizontal
 - * maximal
 - * minimal
 - * 90° to maximum
 - * 90° to minimum
 - * averaged





Particle Gallery



“Looking for a needle in a haystack”





- ✓ Defined filters with *different size-* and *form conditions*
- ✓ *Unlimited number* of particles within the filter conditions are displayed in a *gallery window*
- ✓ *Ability to trace back* to particle movie



Measuring ball bearing beads



✓ Successful search for *fibres* of the tissue used for cleaning of GRADIS






User defined Results

- ★ Considering the *really important fractions* for the evaluation
- ★ Define *new filter conditions* or transfer the filter conditions of the *particle gallery*
- ★ **Re-calculate size** and *shape information* only for particles matching the filter conditions

Examples

- ★ Calculate individual PSDs for the components of a mixture of components that differ in shape
- ★ Calculate the size distribution for perfect spherical particles only
- ★ Eliminate or search for artefacts





Technical Data

Sensor:	QICPIC	1 μm – 20.000 μm
Principle:	Image Analysis	in transmission
Dispersion:	Adaptable modules	for aero-dispersion, sprays, suspensions
Light source:	Pulsed light source	0 to 500 Hz, < 1 ns, 532 nm (green)
	Beam expansion	adaptable to measuring range
Measuring ranges:	Optical modules	3:1 to 1:10, on revolving disk
Measurement:	High speed CMOS camera	1024 x 1024 square pixels, 10 x 10 μm^2 256 grey levels (8 bit)
	Frame rate	up to 500 images/s
Data transmission:	Twisted, two-wire copper cable	1.25 Gbit/s, 5 meters max.
Evaluation:	Particle size	equivalent sphere, enclosing rectangle, cord length, Feret & Martin diameter
	Particle shape	sphericity, as function of particle size
	Class limits	definable, internal 20.000 classes
Visualisation:	Presentation of characteristic particle images	various selection options: movie, gallery
Software & QA:	WINDOX 5	compliant with 21 CFR Rule 11, including electronic signatures



QICPIC Conclusions

- ✓ QICPIC is an innovative measuring system for *particle size* and *shape determination* with *image analysis* in the range of *1 μm* to *20 mm*
- ✓ Extraordinary short exposure time of *< 1 ns* enables clear images of fastest particles and for the first time the use of approved, *cutting edge dispersing devices*, for dry, wet, spray or inhalation applications
- ✓ High acquisition rate of un-matched *500 images/second* guarantees sound particle counts and *superior statistical relevance* of results
- ✓ *Permanent storage* of original information in *20.000 size classes* in data base
- ✓ *Effective data compression* w/o loss of data allows retrospective modifications of evaluation and visualisation
- ✓ Various different *size* and *shape factors* are implemented
- ✓ “*Particle Gallery*” or “*Evaluation of User specific Fractions*” offer powerful tools for investigation of specialities
- ✓ WINDOX 5 software with *client/server structure* and to speed multi-user data base accomplishes *stunning data volume* and manages simultaneous operation with *Laser Diffraction*, *Ultrasonic Extinction* and *Photon Cross Correlation sensors*

