


Quality management, validation of instruments and
global standardisation



Quality

Sympatec's Quality Securing System

SQS


HST: NT 0903-06



Validation of Sympatec Instruments

Introduction


- ★ **Principle** of Laser Diffraction, e.g.
 - ☆ Fraunhofer vs. Mie
- ★ **Verification**
 - ☆ Reference materials vs. reticles
 - ☆ Reference material sets
 - ☆ Specifications
- ★ **Sympatec's Quality Securing System (SQS)**
 - ☆ Staff
 - ☆ Development
 - ☆ Production
 - * Flow
 - * Final control
 - * Certification
 - ☆ Re-certification
 - ☆ Validation binder
 - ☆ Audits



Principle of Laser Diffraction, e.g.


Fraunhofer- vs. Mie-Theory

Physics of interaction



Mie: *Subject to parameters*

Fraunhofer: *Parameter free*
Can be standardised



Validation of LD systems \Rightarrow Fraunhofer



Verification

Reference Material

- ★ *Since 1984* Sympatec applies reference materials for HELOS
- ★ *Verification* of quality of Sympatec instruments
- ★ Market introduction in January 1995
- ★ Backbone of *validation*


Definitions

- ★ **Reference material** = stable material of *arbitrary particle shape*, it is certified values are referenced to a *reference method*
- ★ **Standard material** = *spherical material*, certified values are directly traced to the *standard metre* (via microscopy)

Application Area


- ★ Check of measuring system by comparison with average values of all *Sympatec instruments* of equal configuration
- ★ *Dispersing system* is included in the check






Reference Material Sets

Set	Material	Content	Range
Set 1 HZ0010	SiC-F1200 $x_{50} \approx 4.1 \mu\text{m}$	8 sample bottles of 5 g each 2 bottles of dispersing agent, 10 ml spatula, specifications TKSI-...	R1, R2, R3
Set 2 HZ0020	SiC-P600 $x_{50} \approx 28 \mu\text{m}$	8 sample bottles of 10 g each 2 bottles of dispersing agent, 10 ml spatula, specifications TKSI-...	R3, R4, R5
Set 3 HZ0030	SiC-P50 $x_{50} \approx 420 \mu\text{m}$	8 sample bottles of 20 g each specifications TKSI-...	R6, R7, R8 M7, M8
Set 5 HZ0050	SiC-P80 $x_{50} \approx 250 \mu\text{m}$	8 sample bottles of 20 g each specifications TKSI-...	R5, R6, R7 M6, M7
Set 6 HZ0060	SiC-F1200 $x_{50} \approx 4.1 \mu\text{m}$	64 sample tubes of 1,2 g each, incl. Barcode, specifications TKSI...	ASPIROS R1, R2, R3
Set 7 HZ0070	SiC-P600 $x_{50} \approx 28 \mu\text{m}$	64 sample tubes of 150 mg each, incl. Barcode, specifications TKSI...	ASPIROS R3, R4, R5
Set 8 HZ0080	SiC-P16 $x_{50} \approx 1600 \mu\text{m}$	8 sample bottles of 20 g each specifications TKSI-...	R8 M8




Sympatec's Quality Securing System (SQS)

- ★ Development
 - ☆ Mechanics
 - ☆ Hardware
 - ☆ Software
 - ☆ Firmware
 - ☆ Laser diffraction
 - ☆ Ultrasonic extinction
- ★ Audits
- ★ Staff
- ★ Production
 - ☆ Goods inward
 - ☆ Assembling
 - ☆ Final control
 - ☆ Certification
- ★ Re-certification
- ★ Validation binder



Production



Flow


- ✓ Production file
- ✓ Definition of components
- ✓ Documentation of serial numbers
- ✓ Documented tests
- ✓ Documented responsibility

Sympatec GmbH
System-Produkt-Technik

WINDOX

HELOS (H0841) RODOS+VIBRI: Silikon

VIBRI Measuri
no cascade measur
1,00 bar measur
52,00 mbar cycle
60,00 % start
18,00 % referen
Silicon carbide evaluat
3,22 g/cm3 operato
1,00 identif
KF 813, RO 759, DS 490, Vibri



Volume Size Distribution			
x0/µm	Q3/ %	x0/µm	Q3/ %
0,18	0,02	0,74	3,57
0,22	0,14	0,86	3,83
0,26	0,37	1,00	4,02
0,30	0,66	1,20	4,28
0,36	1,16	1,50	4,99
0,44	1,84	1,80	6,22
0,52	2,46	2,10	8,19
0,62	3,08	2,50	12,31
		8,60	96,28
		10,20	99,23
		35,00	100,00

x10 = 2,28 µm	x50 = 4,45 µm	x90 = 7,36 µm
x16 = 2,74 µm	x84 = 6,88 µm	x99 = 10,07 µm
Sv = 1,9906 m2/cm3	Sm = 6187,7 cm2/g	copt = 4,80 %

Produktionsmappe
ASTRA
H0841/9600 S 220 V / 50 Hz KW 03 / 04 / 06



3 AV Mechanik Dispersiersysteme

RODOS	Bezeichnung	Seriennr.	Fertig	Syst.Int.
<input checked="" type="checkbox"/>	Typ: RODOS 4mm	759	/	/
<input checked="" type="checkbox"/>	Disp.Syst. 4mm gehärtet	DS 490	/	/
<input type="checkbox"/>	Dip.Syst. 4mm nicht geh.			
<input checked="" type="checkbox"/>	Austerhilfe		/	/

für RODOS:				
<input checked="" type="checkbox"/>	Plattform für Tellerbetrieb		/	/
<input checked="" type="checkbox"/>	Teller 1mm	<input type="checkbox"/> Teller 2mm	/	/
<input checked="" type="checkbox"/>	Bürste	(nicht bei 6mm)	/	/
<input checked="" type="checkbox"/>	Ersatzbürste	(nicht bei 6mm)	/	/
<input checked="" type="checkbox"/>	Rüssel		/	/
<input checked="" type="checkbox"/>	Trichteradapter	<input type="checkbox"/> 6mm	/	/
<input checked="" type="checkbox"/>	Plattform für Trichterbetrieb		/	/
<input type="checkbox"/>	Kaskade	(nicht bei 6mm)		
<input type="checkbox"/>	Stromversorgung RODOS			
<input checked="" type="checkbox"/>	RS232C-Kabel		/	/
<input type="checkbox"/>	Vario-Unterbau		/	/
<input type="checkbox"/>	Zusatzzeller 1mm	<input type="checkbox"/> Zusatzzeller 2mm		
<input type="checkbox"/>	Zusatzinjektor	<input type="checkbox"/> 6mm		
<input checked="" type="checkbox"/>	RI-Absaugung		/	/

für auto-RODOS Modul				
<input type="checkbox"/>	Trichteradapter	<input type="checkbox"/> 6mm		
<input type="checkbox"/>	DSD	<input type="checkbox"/> DRB <input type="checkbox"/> DTB		
<input type="checkbox"/>	Vorbeschleunigung			
<input type="checkbox"/>	Trockenmelballe	<input type="checkbox"/> Tactatur (auto-RODOS)		
<input checked="" type="checkbox"/>	Adapter für Absaugung		/	/
<input checked="" type="checkbox"/>	Erdungskabel		/	/
<input checked="" type="checkbox"/>	Schnellverschluss		/	/

GRADIS	Bezeichnung	Seriennr.	Fertig	Syst.Int.
<input type="checkbox"/>	GRADIS	<input type="checkbox"/> Vario		
<input type="checkbox"/>	auto-GRADIS (mit Staubsaugeranschaltbox)			
<input type="checkbox"/>	Stativ			
<input type="checkbox"/>	Plattform für Netzanschluß			
<input type="checkbox"/>	Netzkabel für Plattform			

ZA	Bezeichnung	Seriennr.	Fertig	Syst.Int.
<input type="checkbox"/>	Niifisk + Zyklon+Verschlauchung			
<input type="checkbox"/>	Ersatzstaubbeutel (Niifisk)			
<input type="checkbox"/>	Absolutfilter			
<input type="checkbox"/>	Kärcher			
<input type="checkbox"/>	Verschlauchung + Gummimuffe			
<input type="checkbox"/>	Gewebefilter (Kärcher)			
<input checked="" type="checkbox"/>	Schwingmehrdosierer	191	/	/
<input checked="" type="checkbox"/>	Trichter für Dosierlinie		/	/

Copyright © 1995, Sympatec GmbH, System-Produkt-Technik, Burggrabenstraße 4, D-38678 Clausthal-Zellerfeld, Deutschland

Produktionsmappe
ASTRA
H0841/9600 S 220 V / 50 Hz KW 03 / 04 / 06




1 Kontrollblatt

1. Schubert	Produktionsmappe geneigert
Datum: 06.01.97	Unterschrift: <i>[Signature]</i>
2. Schubert	Arbeitsvorbereitung
Datum: 06.01.97	Unterschrift: <i>[Signature]</i>
3. Schubert	Vorbereitung Software
Datum: 06.01.97	Unterschrift: <i>[Signature]</i>
4. Migge	Vorbereitung Einlauf
Datum: 07.01.97	Unterschrift: <i>[Signature]</i>
5. Isensee	Produktion Mechanik
Datum: 08.01.97	Unterschrift: <i>[Signature]</i>
6. Riedel	Optronik / Hardware
Datum: 09.01.97	Unterschrift: <i>[Signature]</i>
7. Hallenberger	Systemintegration
Datum: 17.01.97	Unterschrift: <i>[Signature]</i>
8. Witt	Freigabe
Datum: 17.01.97	Unterschrift: <i>[Signature]</i>
9. Ludwig	Zusammenfügen mit Auftragsmappe
Datum: 19.01.97	Unterschrift: <i>[Signature]</i>

Copyright © 1995, Sympatec GmbH, System-Produkt-Technik, Burggrabenstraße 4, D-38678 Clausthal-Zellerfeld, Deutschland

Final Control

- ★ Three stages
 - ✓ Test on reception (laser, detector, optics, etc.)
 - ✓ „Cold“ check (mechanical, electrical, optical)
 - ✓ „Hot“ check (measurements with particles)
- ★ Completely documented with production file



- ★ All measuring results available as data-files
- ★ Extent: Currently >> 150.000 files



Certification

Every Sympatec Instrument has to perform a four-step certification procedure during the system integration process:

1. **Electric** and **optical properties** are measured, documented and compared with pre-defined limits
2. Every possible combination of measuring range and dispersing unit is tested with *reference material*:
 - ✓ In general two materials are measured at least three times
 - ✓ X₁₀, X₁₆, X₅₀, X₈₄, X₉₀ (and X₉₉) are listed in the specification sheet
 - ✓ The average values are calculated and compared with the established absolute limits
 - ✓ All serial numbers are taken to be listed in the certification sheet
 - ✓ The production file is completed with all the manufacturing notes, signatures, serial numbers and measured results
3. The certificate is *signed off* only, if all results are within the limits and the documentation is complete
4. The system is shipped together with the certificate



CERTIFICATE

HELOS system no. H1234
comprising:

Component	Serial Number	Component	Serial Number
HELOS/BF	1331	MEGAGRADIS	247
Measuring ranges	R1-R7T	VIBRI	
PCI-FOL	279A	ASPIROS	111
WINDOX 4	Vers. 2 Rel.1.0	SUCCELL	133
MIE	050A	Cuvette 2mm	02025
Add-ons		Cuvette 4mm	02024
USB key 1	2012U	QUIXEL	258
USB key 2		CUV-CHASSIS	
QX		CUV-50ML/US	
HP-Vectra	YBQN140632	Glass 50ml	
		Spare glass 50ml	
WINDOWS 2000		CUV-6ML/SM	
Monitor	YEUN088617	Glass 6ml	
Keyboard	YBKB020617279786	Spare glass 6ml	
		SAFIR	
Mouse	HCA20570877	INCELL	
		SPRAYER	152
Printer		ROTOR/M	
RODOS	1259	INHALER	
Dispersing line 4mm	DS1213	Vacuum pump	
Dispersing line 6mm			
Nilfisk			
Compressor			

The system listed above was carefully manufactured and extensively tested by Sympatec's System Integration Group.

We certify, that the testing results for the preceding configuration meet the Sympatec specifications listed as TKSI02-4DR without any restriction.

Goslar,
SYMPATEC GmbH
System-Partikel-Technik

Goslar,
SYMPATEC GmbH
System-Partikel-Technik

Head of Quality Assurance

Head of System Integration Group

Production Certificate



CERTIFICATE

HELOS system no. H1234
comprising:

Component	Serial Number/Chk	Component	Serial Number/Chk
HELOS/BF	1331	MEGAGRADIS	247
Measuring ranges	R1-R7T	VIBRI	
PCI-FOL	279A	ASPIROS	111
WINDOX 4	Vers. 2 Rel.1.0	SUCCELL	133
MIE	050A	Cuvette 2mm	02025
Add-ons		Cuvette 4mm	02024
USB key 1	2012U	QUIXEL	258
USB key 2		CUV-CHASSIS	
QX		CUV-50ML/US	
HP-Vectra	YBQN140632	Glass 50ml	
		Spare glass 50ml	
WINDOWS 2000		CUV-6ML/SM	
Monitor	YEUN088617	Glass 6ml	
Keyboard	YBKB020617279786	Spare glass 6ml	
		SAFIR	
Mouse	HCA20570877	INCELL	
		SPRAYER	152
Printer		ROTOR/M	
RODOS	1259	INHALER	
Dispersing line 4mm	DS1213	Vacuum pump	
Dispersing line 6mm			
Nilfisk			
Compressor			

We certify, that the testing results for the preceding configuration meet the Sympatec specifications listed as TKSI02-4DR without any restriction.

The system listed above was inspected and tested after installation by the signer


Goslar,
SYMPATEC GmbH
System-Partikel-Technik

Place, date
Company
Signature

Head of Quality Assurance

Name
Inspector

Installation Certificate



Re-certification

Required Number of Reference Material Sets:					
Dispersing Unit	Measuring Range	HZ0010	HZ0020	HZ0030	HZ0050
		SiC-F1200 8 x 10 g	SiC-P600 8 x 10 g	SiC-P50 8 x 20 g	SiC-P80 8 x 20 g
RODOS,	R1 (f = 20 mm)	4			
RODOS/M	R2 (f = 50 mm)	4			
4 mm	R3 (f = 100 mm)	4	4		
Disp.strecke	R4 (f = 200 mm)		4		
	R5 (f = 500 mm)		4		4
	R6 (f = 1000 mm)				4
	R7 (f = 2000 mm)				4
	6 mm Dispersing Line only	R6 (f = 1000 mm) R7 (f = 2000 mm)			4 4
GRADIS, MEGAGRADIS	R5 (f = 500 mm)				4
	R6 (f = 1000 mm)			4	4
	R7 (f = 2000 mm)			4	4
	R8 (f = 5000 mm)			8	
SUCELL, SUCELL/M, QUIXEL	R1 (f = 20 mm)	0.1			
	R2 (f = 50 mm)	0.1			
	R3 (f = 100 mm)	0.1	0.1		
	R4 (f = 200 mm)		0.1		
	R5 (f = 500 mm)		0.1		
	R6 (f = 1000 mm)		0.1		
CUVETTE 50 ml	R2 (f = 50 mm)	0.1			
	R3 (f = 100 mm)	0.1	0.1		
	R4 (f = 200 mm)		0.1		
	R5 (f = 500 mm)		0.1		
	R6 (f = 1000mm)		0.1		

Any Sympatec Instrument can be re-certified at customer's location at any time using reference material sets.

Conditions for approval of re-certification:

- ★ Serial numbers / results in accordance with all limits
- ★ Testing person approves validity of procedure with his signature




Validation Binder

Contents

- ★ General:
 - ☆ Sympatec's validation procedure
 - ☆ Sympatec's SQS strategy
 - ☆ Sample quotation for a validation procedure
 - ☆ Example of an audit procedure
 - ☆ Organisation chart and status of education
 - ☆ Software validation
 - ☆ Back-up strategy
- ★ System specific documents:
 - ☆ Verification of system specifications
 - ☆ Certificate
 - ☆ Copy of production file (in parts)
 - ☆ Specifications of reference materials
- ★ system specific test procedures (IQ, OQ, PQ) for:
 - ☆ Installation
 - ☆ Documentation
 - ☆ Software files
 - ☆ Ambient conditions
 - ☆ Re-certification (prepared)





Conclusion

- ★ Sympatec's validation procedure of instruments for particle size analysis has been established as *standard procedure*
- ★ All systems are produced in accordance to the requests of an *optional subsequent validation*
- ★ The validation procedure is supported with
 - ☆ Sympatec's Validation Binder
 - ☆ Sympatec's Reference Materials
 - ☆ Optional audits of Sympatec's production
 - ☆ Certified service engineers
 - ☆ Adapted service contracts
- ★ Currently validation is available for
 - ☆ HELOS*
 - ☆ MYTOS*
 - ☆ OPUS
 - ☆ QICPIC


Sympatec's *validation procedure* is the standard procedure for installations in pharmaceutical industry

✓ * both are compliant with *DIN/ISO 13320-1*



CFR 21 Part 11

Sympatec's Contribution to Compliance



History

- ★ In 1997, the American Commissioner of Food and Drugs published an amendment to the Code of Federal Regulations, Title 21, Chapter 1, which was entitled “part 11 - Electronic Records; Electronic Signatures”.
- ★ The main objective of what was then called “Rule 11” is to set up regulations for *electronic records* that correspond to the existing regulations for *records on paper*, and thus provide rules for the *paperless laboratory*.
- ★ Since none of the computerised systems available at that time was compliant with this new regulation and all manufacturers claimed time to adapt to the rules, execution of the so called “Rule 11” was suspended for some time.
- ★ Since 2001, the FDA has been urging the pharmaceutical industry to achieve compliance with Rule 11.
 - ↳ Consequently, our pharmaceutical customers urged us for support




What Is This Rule 11 About?

Rule 11 contains regulations for the paperless lab. It takes respect to the trend that more and more information is stored electronically instead of on paper.

It contains

- ★ Rules for the storage, archiving, and restoration of *electronic records*
- ★ Rules to identify the *authorship* and maintain the *integrity* of the information (write-protection of data and/or *audit trail* of modifications)
- ★ Rules how *electronic signatures* can be applied to electronic records and under which circumstances they are considered *equivalent* to a signature on paper





Rule 11 is Optional !

Preamble chapter XVI. *Analysis of Impacts*, Section C reads as follows:

1.1 C. Description of the Impact

For any paper record that an entity is required to keep under existing statutes or FDA regulations, FDA will now accept an electronic record instead of a paper one, as long as the electronic record conforms to the requirements of this rule. FDA will also consider an electronic signature to be equivalent to a handwritten signature if it meets the requirements of this rule.

.....

.....

This action is voluntary; paper records and handwritten signatures are still fully acceptable. No entity will be required to change the way it is currently allowed to submit paper records to the agency.



Where can I find more Information ?

It would exceed the scope of this lecture to present more details. The Rule 11 document is available from the FDA's internet page (PDF document, 38 pages consisting of 35 pages preamble + 3 pages "official" rule)

http://www.fda.gov/ora/compliance_ref/part11/FRs/background/pt11finr.pdf

Outline of the Document

PART 11—ELECTRONIC RECORDS; ELECTRONIC SIGNATURES

Subpart A—General Provisions

- 11.1 Scope
- 11.2 Implementation
- 11.3 Definitions

Subpart B—Electronic Records

- 11.10 Controls for closed systems
- 11.30 Controls for open systems

11.50 Signature manifestations

11.70 Signature/record linking

Subpart C—Electronic Signatures

11.100 General requirements

11.200 Electronic signature components
and controls

11.300 Controls for identification
codes/passwords



How can Compliance be achieved ?

★ A common question is *“The xyz software, is it compliant ?”*

This is a wrong question. It should be:

“How much support for compliance does the computerised system provide?”

What is the difference between both questions?

★ Compliance is not a matter of *Yes or No* but of *better or worse*.

So there is room for competition among suppliers.

★ It is not only the application software (WINDOX in our case) but the *whole computerised system* that must be considered.

★ Not only the *supplier* but also the *pharmaceutical company* (by SOPs and other administrative measures) must contribute to achieve compliance.

↪ **Compliance requests *cooperation* of the *supplier* and the *pharmaceutical company*.**





What is Sympatec's Contribution to Compliance?

★ Sympatec prepared a Rule 11 checklist with:

★ *Answers* to a potential auditor's questions about our systems with respect to compliance

★ *Instructions* how to use the WINDOX Rule 11 features properly

★ *Proposals* of additional measures that the Pharmaceutical company should take

★ About 85 of the points have to be fulfilled by the *supplier*, Sympatec

★ About 115 by the *pharmaceutical company*



WINDOX 4
Electronic Records/Electronic Signatures
Compliance Assessment Worksheet for 21 CFR Part 11

Reference	Question	Assessment Result	Remarks
11.10(b) (continued) Preamble clauses 69-70	• Can a copy of a single record (in electronic format) be supplied to an inspector? In paper format?	S:Yes	Electronic format: Database export function. Paper format: Report.
	• Can a copy of the entire database (in electronic format) be supplied to an inspector?	S:Yes	The database is a directory of files. It can be copied onto any medium with the operating system's file explorer.
	• Are procedures in place to describe HOW to accomplish these inspection tasks?	S:Yes	We recommend to use the Service Request Wizard. It creates a compressed copy of all necessary files.
	• Are procedures in place to define what in format the electronic records will be provided?	S:Yes	On-line-Help contains a description of all Report template statements, their value and format.
11.10(c) Preamble clause 71	• Are the records protected to ensure their accurate and ready retrieval throughout the record retention period?	S:Yes	No measured data or data describing the circumstances of the measurement can be deleted from the database. If the database format is changed due to an update, a conversion utility is provided to ensure compatibility of old data.
	• Are records protected on the system to prevent unauthorized modification or	S:Yes	



How does the WINDOX Software support Rule 11 ?

★ System requirements:

To ensure the correctness of the user identification, an *operating system with access control* (Windows NT, 2000, or XP) must be used.

★ Standard features:

These features are supplied with *WINDOX* and are always enabled.

- ☆ WINDOX has been developed under conditions of *Life Cycle Development Quality Control*. This is a must for Rule 11 compliance.
- ☆ WINDOX stores *user logon, date and time* with every data record.
- ☆ WINDOX keeps an *administration log* in which all administrative action (e.g. disabling of "Rule 11 mode") is written.
- ☆ *Data integrity* is enforced by the WINDOX database.
- ☆ WINDOX is *downward compatible* to the respective predecessor generation of data records. This ensures the restoration of archived data, even if, sometimes, in consecutive steps (HELOS/DOS → WINDOX 3 → WINDOX 4 → WINDOX 5).




How does the WINDOX Software support Rule 11 ?

(continued)

★ Additional measures:


These features are supplied but disabled by default, or options to call.

☆ WINDOX has a “*Rule 11 mode*” of database operation and several more switches to enforce data integrity.



☆ WINDOX comes with an optional *digital signature*.

Measurements can be reviewed and signed by someone other than the operator.



Conclusions

Current State

↪ *WINDOX fulfils the supplier's part of Rule 11 compliance with respect to **digital records**, has improved features of data protection and audit trail recording and comes with a **digital signature***

Future Steps

- ★ Sympatec will now wait for *feedback* from the Pharmaceutical Industry about the benefits and eventual problems of our Rule 11 support
- ★ A new challenge is the proper *archiving and restoration* of electronic records over a retention period of ten years, as required by the FDA

