



**Neutron<sup>®</sup> Pharmacchemical Co.**  
Manufacturer of Laboratory Chemical & Pharmaceutical Materials

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Revision Date 14.07.2018

Version 8.9

---

## SECTION 1. Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

REACH Registration Number This product is a mixture. REACH Registration Number see section 3.

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
	For additional information on uses please refer to the Neutron Pharma Chemicals portal ( <a href="http://www.Neutronpharmchemical.com">www.Neutronpharmchemical.com</a> ).

### 1.3 Details of the supplier of the safety data sheet

Company	NEUTRON PHARMACHEMICAL CO 98, 9th Floor, Borjsaz Building, Azadi Ave, Tehran, Iran.
Responsible Department	T 021-66906732-3 - F 021-66581408 - <a href="mailto:info@neutronco.com">info@neutronco.com</a> <a href="http://www.neutronco.com">www.neutronco.com</a>

1.4 Emergency telephone number	125
--------------------------------	-----

---

## SECTION 2. Hazards identification

### 2.1 Classification of the substance or mixture

Classification (REGULATION (EC) No 1272/2008)

Flammable liquid, Category 2, H225

For the full text of the H-Statements mentioned in this Section, see Section 16.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

## 2.2 Label elements

### Labelling (REGULATION (EC) No 1272/2008)

#### *Hazard pictograms*



#### *Signal word*

Danger

#### *Hazard statements*

H225 Highly flammable liquid and vapour.

#### *Precautionary statements*

##### Prevention

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P240 Ground/bond container and receiving equipment.

##### Storage

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

### Reduced labelling ( $\leq 125$ ml)

#### *Hazard pictograms*



#### *Signal word*

Danger

#### *Precautionary statements*

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

## 2.3 Other hazards

None known.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760  
Product name Mixed indicator 4.5 acc. to Mortimer

---

---

## SECTION 3. Composition/information on ingredients

Chemical nature Aqueous-ethanolic dye solution.

### 3.1 Substance

Not applicable

### 3.2 Mixture

#### Hazardous components (REGULATION (EC) No 1272/2008)

*Chemical name (Concentration)*

CAS-No.	Registration number	Classification
---------	---------------------	----------------

ethanol ( $\geq 25\%$ - $< 50\%$ )		
------------------------------------	--	--

*Substance does not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII.*

64-17-5	01-2119457610-43-XXXX	Flammable liquid, Category 2, H225 Eye irritation, Category 2, H319
---------	-----------------------	--

For the full text of the H-Statements mentioned in this Section, see Section 16.

---

## SECTION 4. First aid measures

### 4.1 Description of first aid measures

After inhalation: fresh air.

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

After eye contact: rinse out with plenty of water.

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

### 4.2 Most important symptoms and effects, both acute and delayed

irritant effects, Nausea, Vomiting, inebriation, euphoria, respiratory paralysis, Dermatitis, Dizziness, narcosis

### 4.3 Indication of any immediate medical attention and special treatment needed

No information available.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

---

## SECTION 5. Firefighting measures

### 5.1 Extinguishing media

#### *Suitable extinguishing media*

Water, Foam, Carbon dioxide (CO<sub>2</sub>), Dry powder

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### 5.2 Special hazards arising from the substance or mixture

Combustible.

Vapours are heavier than air and may spread along floors.

Pay attention to flashback.

Forms explosive mixtures with air at ambient temperatures.

Development of hazardous combustion gases or vapours possible in the event of fire.

### 5.3 Advice for firefighters

#### *Special protective equipment for firefighters*

In the event of fire, wear self-contained breathing apparatus.

#### *Further information*

Remove container from danger zone and cool with water. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## SECTION 6. Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Do not breathe vapours, aerosols. Ensure adequate ventilation. Keep away from heat and sources of ignition. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders:

Protective equipment see section 8.

### 6.2 Environmental precautions

Do not let product enter drains. Risk of explosion.

### 6.3 Methods and materials for containment and cleaning up

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760  
Product name Mixed indicator 4.5 acc. to Mortimer

---

Cover drains. Collect, bind, and pump off spills. Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent and neutralising material (e.g. Chemizorb® H<sup>+</sup>, Merck Art. No. 101595). Dispose of properly. Clean up affected area.

## 6.4 Reference to other sections

Indications about waste treatment see section 13.

---

## SECTION 7. Handling and storage

### 7.1 Precautions for safe handling

#### *Advice on safe handling*

Observe label precautions.

#### *Advice on protection against fire and explosion*

Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharge.

#### *Hygiene measures*

Change contaminated clothing. Preventive skin protection recommended. Wash hands after working with substance.

### 7.2 Conditions for safe storage, including any incompatibilities

#### *Storage conditions*

Keep container tightly closed in a dry and well-ventilated place. Keep away from heat and sources of ignition.

Recommended storage temperature see product label.

### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

---

## SECTION 8. Exposure controls/personal protection

### 8.1 Control parameters

#### Derived No Effect Level (DNEL)

##### *ethanol (64-17-5)*

Worker DNEL, acute	Local effects	inhalation	1900 mg/m <sup>3</sup>
--------------------	---------------	------------	------------------------

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760  
Product name Mixed indicator 4.5 acc. to Mortimer

---

Worker DNEL, longterm	Systemic effects	dermal	343 mg/kg Body weight
Worker DNEL, longterm	Systemic effects	inhalation	950 mg/m <sup>3</sup>
Consumer DNEL, acute	Local effects	inhalation	950 mg/m <sup>3</sup>
Consumer DNEL, longterm	Systemic effects	dermal	206 mg/kg Body weight
Consumer DNEL, longterm	Systemic effects	inhalation	114 mg/m <sup>3</sup>
Consumer DNEL, longterm	Systemic effects	oral	87 mg/kg Body weight

## Predicted No Effect Concentration (PNEC)

*ethanol (64-17-5)*

PNEC Fresh water	0,96 mg/l
PNEC Marine water	0,79 mg/l
PNEC Fresh water sediment	3,6 mg/kg
PNEC Soil	0,63 mg/kg
PNEC Aquatic intermittent release	2,75 mg/l
PNEC Sewage treatment plant	580 mg/l
PNEC oral	720 mg/kg

## 8.2 Exposure controls

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

See section 7.1.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760  
Product name Mixed indicator 4.5 acc. to Mortimer

---

## Individual protection measures

Protective clothing needs to be selected specifically for the workplace, depending on concentrations and quantities of the hazardous substances handled. The chemical resistance of the protective equipment should be enquired at the respective supplier.

### *Eye/face protection*

Safety glasses

### *Hand protection*

full contact:

Glove material:	butyl-rubber
Glove thickness:	0,7 mm
Break through time:	> 480 min

splash contact:

Glove material:	Nitrile rubber
Glove thickness:	0,40 mm
Break through time:	> 120 min

The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374, for example KCL 898 Butoject® (full contact), KCL 730 Camatril® -Velours (splash contact).

The breakthrough times stated above were determined by KCL in laboratory tests acc. to EN374 with samples of the recommended glove types.

This recommendation applies only to the product stated in the safety data sheet(>,<)> supplied by us and for the designated use. When dissolving in or mixing with other substances and under conditions deviating from those stated in EN374 please contact the supplier of CE-approved gloves (e.g. KCL GmbH, D-36124 Eichenzell, Internet: [www.kcl.de](http://www.kcl.de)).

### *Other protective equipment*

Flame retardant antistatic protective clothing.

### *Respiratory protection*

required when vapours/aerosols are generated.

Recommended Filter type: filter ABEK

The entrepreneur has to ensure that maintenance, cleaning and testing of respiratory protective devices are carried out according to the instructions of the producer. These measures have to be properly documented.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

## Environmental exposure controls

Do not let product enter drains.

Risk of explosion.

---

## SECTION 9. Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

Form	liquid
Colour	red brown
Odour	of ethanol
Odour Threshold	No information available.
pH	ca. 4,8 at 20 °C
Melting point	No information available.
Boiling point	No information available.
Flash point	9 °C
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapour pressure	No information available.



# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

Relative vapour density	No information available.
-------------------------	---------------------------

Density	0,93 g/cm <sup>3</sup> at 20 °C
---------	------------------------------------

Relative density	No information available.
------------------	---------------------------

Water solubility	at 20 °C soluble
------------------	---------------------

Partition coefficient: n-octanol/water	No information available.
--	---------------------------

Auto-ignition temperature	No information available.
---------------------------	---------------------------

Decomposition temperature	No information available.
---------------------------	---------------------------

Viscosity, dynamic	No information available.
--------------------	---------------------------

Explosive properties	Not classified as explosive.
----------------------	------------------------------

Oxidizing properties	none
----------------------	------

## 9.2 Other data

none

---

## SECTION 10. Stability and reactivity

### 10.1 Reactivity

Vapours may form explosive mixture with air.

### 10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### 10.3 Possibility of hazardous reactions

Risk of explosion/exothermic reaction with:

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

hydrogen peroxide, perchlorates, perchloric acid, Nitric acid, mercury(II) nitrate, permanganic acid, Nitriles, peroxi compounds, Strong oxidizing agents, nitrosyl compounds, Peroxides, sodium, Potassium, halogen oxides, calcium hypochlorite, nitrogen dioxide, metallic oxides, uranium hexafluoride, iodides, Chlorine, Alkali metals, Alkaline earth metals, alkali oxides, Ethylene oxide

silver, with, Nitric acid

silver compounds, with, Ammonia

potassium permanganate, with, conc. sulfuric acid

Risk of ignition or formation of inflammable gases or vapours with:

halogen-halogen compounds, chromium(VI) oxide, chromyl chloride, Fluorine, hydrides, Oxides of phosphorus, platinum

Nitric acid, with, potassium permanganate

## 10.4 Conditions to avoid

Warming.

## 10.5 Incompatible materials

rubber, various plastics

## 10.6 Hazardous decomposition products

no information available

---

## SECTION 11. Toxicological information

### 11.1 Information on toxicological effects

#### Mixture

*Acute oral toxicity*

Symptoms: Nausea, Vomiting

*Acute inhalation toxicity*

Symptoms: slight mucosal irritations

absorption

*Acute dermal toxicity*

This information is not available.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

## *Skin irritation*

After long-term exposure to the chemical: Dermatitis

## *Eye irritation*

This information is not available.

## *Sensitisation*

This information is not available.

## *Germ cell mutagenicity*

This information is not available.

## *Carcinogenicity*

This information is not available.

## *Reproductive toxicity*

This information is not available.

## *Teratogenicity*

This information is not available.

## *Specific target organ toxicity - single exposure*

This information is not available.

## *Specific target organ toxicity - repeated exposure*

This information is not available.

## *Aspiration hazard*

This information is not available.

## **11.2 Further information**

Systemic effects:

euphoria

After absorption of large quantities:

Dizziness, inebriation, narcosis, respiratory paralysis

Other dangerous properties can not be excluded.

Handle in accordance with good industrial hygiene and safety practice.

## **Components**

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

## *ethanol*

### *Acute oral toxicity*

LD50 Rat: 10.470 mg/kg

OECD Test Guideline 401

### *Acute inhalation toxicity*

LC50 Rat: 124,7 mg/l; 4 h ; vapour

OECD Test Guideline 403

### *Skin irritation*

Rabbit

Result: No skin irritation

OECD Test Guideline 404

### *Eye irritation*

Rabbit

Result: Eye irritation

OECD Test Guideline 405

### *Sensitisation*

Sensitisation test (Magnusson and Kligman):

Result: negative

(IUCLID)

### *Germ cell mutagenicity*

#### *Genotoxicity in vitro*

Ames test

Salmonella typhimurium

Result: negative

Method: OECD Test Guideline 471

In vitro mammalian cell gene mutation test

Mouse lymphoma test

Result: negative

Method: OECD Test Guideline 476

### *Reproductive toxicity*

Application Route: Oral

Mouse

Method: OECD Test Guideline 416

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

---

## SECTION 12. Ecological information

### Mixture

#### 12.1 Toxicity

No information available.

#### 12.2 Persistence and degradability

No information available.

#### 12.3 Bioaccumulative potential

No information available.

#### 12.4 Mobility in soil

No information available.

#### 12.5 Results of PBT and vPvB assessment

Substance(s) in the mixture do(es) not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII, or a PBT/vPvB assessment was not conducted.

#### 12.6 Other adverse effects

Discharge into the environment must be avoided.

### Components

#### *ethanol*

##### *Toxicity to fish*

LC50 *Leuciscus idus* (Golden orfe): 8.140 mg/l; 48 h  
(IUCLID)

##### *Toxicity to daphnia and other aquatic invertebrates*

EC5 *E.sulcatum*: 65 mg/l; 72 h  
(Lit.)

EC50 *Daphnia magna* (Water flea): 9.268 - 14.221 mg/l; 48 h  
(IUCLID)

##### *Toxicity to algae*

IC5 *Scenedesmus quadricauda* (Green algae): 5.000 mg/l; 7 d  
(Lit.)

##### *Toxicity to bacteria*

EC5 *Pseudomonas putida*: 6.500 mg/l; 16 h  
(IUCLID)

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

*Toxicity to daphnia and other aquatic invertebrates (Chronic toxicity)*

semi-static test NOEC Daphnia magna (Water flea): 9,6 mg/l; 9 d

(ECHA)

*Biodegradability*

94 %

OECD Test Guideline 301E

Readily biodegradable

*Biochemical Oxygen Demand (BOD)*

930 - 1.670 mg/g (5 d)

(Lit.)

*Theoretical oxygen demand (ThOD)*

2.100 mg/g

(Lit.)

*Ratio COD/ThBOD*

90 %

(Lit.)

*Partition coefficient: n-octanol/water*

log Pow: -0,31

(experimental)

(Lit.) Bioaccumulation is not expected.

Substance does not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

---

## SECTION 13. Disposal considerations

### *Waste treatment methods*

See [www.retrologistik.com](http://www.retrologistik.com) for processes regarding the return of chemicals and containers, or contact us there if you have further questions.

---

## SECTION 14. Transport information

### Land transport (ADR/RID)

14.1 UN number	UN 1170
14.2 Proper shipping name	ETHANOL
14.3 Class	3
14.4 Packing group	II
14.5 Environmentally hazardous	--
14.6 Special precautions for user	yes
Tunnel restriction code	D/E

### Inland waterway transport (ADN)

Not relevant

### Air transport (IATA)

14.1 UN number	UN 1170
14.2 Proper shipping name	ETHANOL
14.3 Class	3
14.4 Packing group	II
14.5 Environmentally hazardous	--
14.6 Special precautions for user	no

### Sea transport (IMDG)

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

14.1 UN number	UN 1170
14.2 Proper shipping name	ETHANOL
14.3 Class	3
14.4 Packing group	II
14.5 Environmentally hazardous	--
14.6 Special precautions for user	yes
EmS	F-E S-D

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code  
Not relevant

---

## SECTION 15. Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### *EU regulations*

Major Accident Hazard	SEVESO III
Legislation	FLAMMABLE LIQUIDS P5c Quantity 1: 5.000 t Quantity 2: 50.000 t

Occupational restrictions	Take note of Dir 94/33/EC on the protection of young people at work.
---------------------------	--

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer	not regulated
---	---------------

Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC	not regulated
---	---------------



# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 1.3760

Product name Mixed indicator 4.5 acc. to Mortimer

---

Substances of very high concern (SVHC)

This product does not contain substances of very high concern according to Regulation (EC) No 1907/2006 (REACH), Article 57 above the respective regulatory concentration limit of  $\geq 0.1\%$  (w/w).

*National legislation*

Storage class 3

## 15.2 Chemical safety assessment

For this product a chemical safety assessment was not carried out.

---

## SECTION 16. Other information

**Full text of H-Statements referred to under sections 2 and 3.**

H225 Highly flammable liquid and vapour.

### Training advice

Provide adequate information, instruction and training for operators.

### Labelling

*Hazard pictograms*



*Signal word*

Danger

*Hazard statements*

H225 Highly flammable liquid and vapour.

*Precautionary statements*

Prevention

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	1.3760
Product name	Mixed indicator 4.5 acc. to Mortimer

---

P240 Ground/bond container and receiving equipment.

Storage

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

## Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

## Regional representation

This information is given on the authorised Safety Data Sheet for your country.

---

*The information contained herein is based on the present state of our knowledge. It characterises the product with regard to the appropriate safety precautions. It does not represent a guarantee of any properties of the product.*