



## SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

### 1.1. Product identifier

|                           |  |
|---------------------------|--|
| Product Description:      | <u>Chromium(III) chloride, anhydrous</u> |
| Cat No. :                 | 12336                                    |
| Synonyms                  | Chromic chloride                         |
| CAS No                    | 10025-73-7                               |
| EC No                     | 233-038-3                                |
| Molecular Formula         | Cl <sub>3</sub> Cr                       |
| REACH registration number | -  |

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

|                      |                          |
|----------------------|--------------------------|
| Recommended Use      | Laboratory chemicals.    |
| Uses advised against | No Information available |

### 1.3. Details of the supplier of the safety data sheet

|         |  |
|---------|--|
| Company | Avocado Research Chemicals Ltd.<br>(Part of Thermo Fisher Scientific)<br>Shore Road, Heysham<br>Lancashire, LA3 2XY,<br>United Kingdom<br>Office Tel: +44 (0) 1524 850506<br>Office Fax: +44 (0) 1524 850608 |
|---------|--|

|                |                                |
|----------------|--------------------------------|
| E-mail address | begel.sdsdesk@thermofisher.com |
|----------------|--------------------------------|

### 1.4. Emergency telephone number

For information **US** call: 001-800-227-6701 / **Europe** call: +32 14 57 52 11  
Emergency Number **US**:001-201-796-7100 / **Europe**: +32 14 57 52 99  
**CHEMTREC** Tel. No. **US**:001-800-424-9300 / **Europe**:001-703-527-3887

## SECTION 2: HAZARDS IDENTIFICATION

### 2.1. Classification of the substance or mixture

CLP Classification - According to GB-CLP Regulations UK SI 2019/720 and UK SI 2020/1567

#### Physical hazards

Substances/mixtures corrosive to metal

Category 1 (H290)

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

## Health hazards

Acute oral toxicity  
Skin Sensitization

Category 4 (H302)  
Category 1 (H317)

## Environmental hazards

Chronic aquatic toxicity

Category 2 (H411)

Full text of Hazard Statements: see section 16

## 2.2. Label elements



Signal Word

Warning

## Hazard Statements

- H290 - May be corrosive to metals
- H302 - Harmful if swallowed
- H317 - May cause an allergic skin reaction
- H411 - Toxic to aquatic life with long lasting effects

## Precautionary Statements

- P390 - Absorb spillage to prevent material damage
- P301 + P330 + P331 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting
- P312 - Call a POISON CENTER or doctor if you feel unwell
- P280 - Wear protective gloves/protective clothing
- P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
- P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention

## 2.3. Other hazards

In accordance with Annex XIII of the REACH Regulation, inorganic substances do not require assessment

Toxic to terrestrial vertebrates

This product does not contain any known or suspected endocrine disruptors

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1. Substances

| Component        | CAS No     | EC No             | Weight % | CLP Classification - According to GB-CLP Regulations UK SI 2019/720 and UK SI 2020/1567       |
|------------------|------------|-------------------|----------|---|
| Chromic chloride | 10025-73-7 | EEC No. 233-038-3 | >95      | Acute Tox. 4 (H302)<br>Skin Sens. 1 (H317)<br>Met. Corr. 1 (H290)<br>Aquatic Chronic 2 (H411) |

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

|                           |   |
|---------------------------|---|
| REACH registration number | - |
|---------------------------|---|

Full text of Hazard Statements: see section 16

## SECTION 4: FIRST AID MEASURES

### 4.1. Description of first aid measures

|                                    |  |
|------------------------------------|--|
| General Advice                     | If symptoms persist, call a  |
| Eye Contact                        | Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.                                  |
| Skin Contact                       | Wash off immediately with plenty of water for at least 15 minutes. If skin irritation persists, call a physician.                                |
| Ingestion                          | Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.  |
| Inhalation                         | Remove to fresh air. If not breathing, give artificial respiration. Get medical attention if symptoms occur.                                     |
| Self-Protection of the First Aider | Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. |

### 4.2. Most important symptoms and effects, both acute and delayed

None reasonably foreseeable. May cause allergic skin reaction. Ingestion causes severe swelling, severe damage to the delicate tissue and danger of perforation: Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain or flushing

### 4.3. Indication of any immediate medical attention and special treatment needed

Notes to Physician                      Treat symptomatically.

## SECTION 5: FIREFIGHTING MEASURES

### 5.1. Extinguishing media

#### Suitable Extinguishing Media

Carbon dioxide (CO<sub>2</sub>). Dry chemical. Chemical foam.

**Extinguishing media which must not be used for safety reasons** No information available.

### 5.2. Special hazards arising from the substance or mixture

Thermal decomposition can lead to release of irritating gases and vapors.

#### Hazardous Combustion Products

Hydrogen chloride gas.

### 5.3. Advice for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## SECTION 6: ACCIDENTAL RELEASE MEASURES

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

## 6.1. Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Use personal protective equipment as required. Avoid dust formation.

## 6.2. Environmental precautions

Do not flush into surface water or sanitary sewer system.

## 6.3. Methods and material for containment and cleaning up

Sweep up and shovel into suitable containers for disposal. Keep in suitable, closed containers for disposal.

## 6.4. Reference to other sections

Refer to protective measures listed in Sections 8 and 13.

## SECTION 7: HANDLING AND STORAGE

### 7.1. Precautions for safe handling

Do not get in eyes, on skin, or on clothing. Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid ingestion and inhalation. Avoid dust formation.

### **Hygiene Measures**

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Wash hands before breaks and after work.

### 7.2. Conditions for safe storage, including any incompatibilities

Corrosives area. Keep containers tightly closed in a dry, cool and well-ventilated place. Store under an inert atmosphere. Protect from moisture.

**Technical Rules for Hazardous Substances (TRGS)** Class 13  
**510 Storage Class (LGK) (Germany)**

### 7.3. Specific end use(s)

Use in laboratories

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1. Control parameters

#### **Exposure limits**

List source(s): **UK** - EH40/2005 Work Exposure Limits, Fourth edition. Published 2020.

| Component        | The United Kingdom  | European Union | Ireland |
|------------------|---|----------------|---------|
| Chromic chloride | STEL: 1.5 mg/m <sup>3</sup> 15 min<br>TWA: 0.5 mg/m <sup>3</sup> 8 hr |                |         |

#### **Biological limit values**

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

#### **Derived No Effect Level (DNEL) / Derived Minimum Effect Level (DMEL)**

No information available

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

## Predicted No Effect Concentration (PNEC)

No information available.

## 8.2. Exposure controls

### Engineering Measures

Ensure adequate ventilation, especially in confined areas. Ensure that eyewash stations and safety showers are close to the workstation location.

Wherever possible, engineering control measures such as the isolation or enclosure of the process, the introduction of process or equipment changes to minimise release or contact, and the use of properly designed ventilation systems, should be adopted to control hazardous materials at source

### Personal protective equipment

#### Eye Protection

Goggles (European standard - EN 166)

#### Hand Protection

Protective gloves

| Glove material | Breakthrough time | Glove thickness | EU standard | Glove comments        |
|----------------|-------------------|-----------------|-------------|-----------------------|
| Natural rubber | See manufacturers | -               | EN 374      | (minimum requirement) |
| Nitrile rubber | recommendations   |                 |             |                       |
| Neoprene       |                   |                 |             |                       |
| PVC            |                   |                 |             |                       |

#### Skin and body protection

Long sleeved clothing.

Inspect gloves before use.

Please observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information)

Ensure gloves are suitable for the task: Chemical compatibility, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion.

Remove gloves with care avoiding skin contamination.

#### Respiratory Protection

When workers are facing concentrations above the exposure limit they must use appropriate certified respirators.

To protect the wearer, respiratory protective equipment must be the correct fit and be used and maintained properly

#### Large scale/emergency use

Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced **Recommended Filter type:** Particulates filter conforming to EN 143

#### Small scale/Laboratory use

Use a NIOSH/MSHA or European Standard EN 149:2001 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.

**Recommended half mask:-** Particle filtering: EN149:2001 When RPE is used a face piece Fit Test should be conducted

#### Environmental exposure controls

Prevent product from entering drains. Do not allow material to contaminate ground water system.

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

### 9.1. Information on basic physical and chemical properties

#### Physical State

Powder Solid

#### Appearance

Purple

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

|  |                          |  |
|--|--------------------------|--|
| <b>Odor</b>                                    | Odorless                 |  |
| <b>Odor Threshold</b>                          | No data available        |  |
| <b>Melting Point/Range</b>                     | 1152 °C / 2105.6 °F      |  |
| <b>Softening Point</b>                         | No data available        |  |
| <b>Boiling Point/Range</b>                     | No information available |  |
| <b>Flammability (liquid)</b>                   | Not applicable           | Solid                                    |
| <b>Flammability (solid,gas)</b>                | No information available |  |
| <b>Explosion Limits</b>                        | No data available        |  |
| <b>Flash Point</b>                             | No information available | <b>Method</b> - No information available |
| <b>Autoignition Temperature</b>                | No data available        |  |
| <b>Decomposition Temperature</b>               | No data available        |  |
| <b>pH</b>                                      | No information available | 2.4 @ 20°C (0.2M)                        |
| <b>Viscosity</b>                               | Not applicable           | Solid                                    |
| <b>Water Solubility</b>                        | Soluble 585 g/L @ 25 °C  |  |
| <b>Solubility in other solvents</b>            | No information available |  |
| <b>Partition Coefficient (n-octanol/water)</b> |                          |  |
| <b>Component</b>                               | <b>log Pow</b>           |  |
| Chromic chloride                               | -3                       |  |
| <b>Vapor Pressure</b>                          | No information available |  |
| <b>Density / Specific Gravity</b>              | 2.8                      |  |
| <b>Bulk Density</b>                            | No data available        |  |
| <b>Vapor Density</b>                           | Not applicable           | Solid                                    |
| <b>Particle characteristics</b>                | No data available        |  |

## 9.2. Other information

|                          |                        |
|--------------------------|------------------------|
| <b>Molecular Formula</b> | Cl <sub>3</sub> Cr     |
| <b>Molecular Weight</b>  | 158.36                 |
| <b>Evaporation Rate</b>  | Not applicable - Solid |

## SECTION 10: STABILITY AND REACTIVITY

### 10.1. Reactivity

None known, based on information available

### 10.2. Chemical stability

Hygroscopic.

### 10.3. Possibility of hazardous reactions

|                                 |                               |
|---------------------------------|-------------------------------|
| <b>Hazardous Polymerization</b> | No information available.     |
| <b>Hazardous Reactions</b>      | None under normal processing. |

### 10.4. Conditions to avoid

Incompatible products. Exposure to moisture. Exposure to moist air or water.

### 10.5. Incompatible materials

Strong oxidizing agents.

### 10.6. Hazardous decomposition products

Hydrogen chloride gas.

## SECTION 11: TOXICOLOGICAL INFORMATION

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Product Information

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

**(a) acute toxicity;**

**Oral** Category 4  
**Dermal** Based on available data, the classification criteria are not met  
**Inhalation** Based on available data, the classification criteria are not met

| Component        | LD50 Oral              | LD50 Dermal             | LC50 Inhalation                    |
|------------------|------------------------|-------------------------|------------------------------------|
| Chromic chloride | LD50 = 440 mg/kg (Rat) | LD50 > 2000 mg/kg (Rat) | 31.5 mg/m <sup>3</sup> /2h (Mouse) |

**(b) skin corrosion/irritation;**

**Test method** OECD 404  
**Test species** rabbit  
**Observational endpoint** No skin irritation

**(c) serious eye damage/irritation;**

Based on available data, the classification criteria are not met  
**Test method** OECD 405  
**Test species** rabbit  
**Observation end point** No eye irritation

**(d) respiratory or skin sensitization;**

**Respiratory** No data available  
**Skin** Category 1

| Component                              | Test method                        | Test species | Study result  |
|--|------------------------------------|--------------|---------------|
| Chromic chloride<br>10025-73-7 ( >95 ) | in vivo<br>OECD Test Guideline 406 | guinea pig   | Sensitization |

No information available

**(e) germ cell mutagenicity;**

Based on available data, the classification criteria are not met

| Component                              | Test method             | Test species | Study result |
|--|-------------------------|--------------|--------------|
| Chromic chloride<br>10025-73-7 ( >95 ) | OECD Test Guideline 473 | in vitro     | negative     |

**(f) carcinogenicity;**

Based on available data, the classification criteria are not met

| Component                              | Test method | Test species / Duration | Study result |
|--|-------------|-------------------------|--------------|
| Chromic chloride<br>10025-73-7 ( >95 ) | in vivo     | Rat                     | negative     |

There are no known carcinogenic chemicals in this product

**(g) reproductive toxicity;**

Based on available data, the classification criteria are not met

| Component                              | Test method             | Test species / Duration | Study result |
|--|-------------------------|-------------------------|--------------|
| Chromic chloride<br>10025-73-7 ( >95 ) | OECD Test Guideline 414 | mouse<br>17 days        | negative     |

**(h) STOT-single exposure;**

No data available

**(i) STOT-repeated exposure;**

No data available

**Target Organs** None known.

**(j) aspiration hazard;**

Not applicable  
Solid

**Other Adverse Effects**

The toxicological properties have not been fully investigated.

**Symptoms / effects, both acute and** Ingestion causes severe swelling, severe damage to the delicate tissue and danger of

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

delayed

perforation. Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain or flushing.

## 11.2. Information on other hazards

### Endocrine Disrupting Properties

Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

## SECTION 12: ECOLOGICAL INFORMATION

### 12.1. Toxicity

#### Ecotoxicity effects

Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment. The product contains following substances which are hazardous for the environment.

| Component        | Freshwater Fish   | Water Flea                              | Freshwater Algae                                 |
|------------------|---|---|--|
| Chromic chloride | LD50 = 57.4 mg/L (96h)<br>Rainbow trout<br>EC10 = 0.246 mg/L<br>Salmo gairdneri | LC50 = 63.3 mg/L (48h)<br>Daphnia magna | EC50 = 2 mg/L (96h)<br>Selenastrum capricornutum |

| Component        | Microtox        | M-Factor |
|------------------|-----------------|----------|
| Chromic chloride | EC50 = 256 mg/L |          |

### 12.2. Persistence and degradability

#### Persistence

May persist, based on information available.

#### Degradability

Not relevant for inorganic substances.

#### Degradation in sewage treatment plant

Contains substances known to be hazardous to the environment or not degradable in waste water treatment plants.

### 12.3. Bioaccumulative potential

May have some potential to bioaccumulate

| Component        | log Pow | Bioconcentration factor (BCF) |
|------------------|---------|-------------------------------|
| Chromic chloride | -3      | No data available             |

### 12.4. Mobility in soil

The product is water soluble, and may spread in water systems Will likely be mobile in the environment due to its water solubility. Highly mobile in soils

### 12.5. Results of PBT and vPvB assessment

In accordance with Annex XIII of the REACH Regulation, inorganic substances do not require assessment.

### 12.6. Endocrine disrupting properties

#### Endocrine Disruptor Information

This product does not contain any known or suspected endocrine disruptors

### 12.7. Other adverse effects

#### Persistent Organic Pollutant

This product does not contain any known or suspected substance

#### Ozone Depletion Potential

This product does not contain any known or suspected substance

## SECTION 13: DISPOSAL CONSIDERATIONS

### 13.1. Waste treatment methods

#### Waste from Residues/Unused Products

Waste is classified as hazardous. Dispose of in accordance with the European Directives on waste and hazardous waste. Dispose of in accordance with local regulations.



# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

**Contaminated Packaging** Dispose of this container to hazardous or special waste collection point.

**European Waste Catalogue (EWC)** According to the European Waste Catalog, Waste Codes are not product specific, but application specific.

**Other Information** Do not flush to sewer. Waste codes should be assigned by the user based on the application for which the product was used. Do not empty into drains. Do not let this chemical enter the environment.

## SECTION 14: TRANSPORT INFORMATION

### IMDG/IMO

**14.1. UN number** UN3260  
**14.2. UN proper shipping name** Corrosive solid, acidic, inorganic, n.o.s.  
**Technical Shipping Name** Chromium(III) chloride  
**14.3. Transport hazard class(es)** 8  
**14.4. Packing group** III

### ADR

**14.1. UN number** UN3260  
**14.2. UN proper shipping name** Corrosive solid, acidic, inorganic, n.o.s.  
**Technical Shipping Name** Chromium(III) chloride  
**14.3. Transport hazard class(es)** 8  
**14.4. Packing group** III

### IATA

**14.1. UN number** UN3260  
**14.2. UN proper shipping name** Corrosive solid, acidic, inorganic, n.o.s.  
**Technical Shipping Name** Chromium(III) chloride  
**14.3. Transport hazard class(es)** 8  
**14.4. Packing group** III

**14.5. Environmental hazards** Dangerous for the environment  
Product is a marine pollutant according to the criteria set by IMDG/IMO

**14.6. Special precautions for user** No special precautions required.

**14.7. Maritime transport in bulk according to IMO instruments** Not applicable, packaged goods

## SECTION 15: REGULATORY INFORMATION

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### International Inventories

Europe (EINECS/ELINCS/NLP), China (IECSC), Taiwan (TCSI), Korea (KECL), Japan (ENCS), Japan (ISHL), Canada (DSL/NDSL), Australia (AICS), New Zealand (NZIoC), Philippines (PICCS). US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

| Component        | CAS No     | EINECS    | ELINCS | NLP | IECSC | TCSI | KECL     | ENCS | ISHL |
|------------------|------------|-----------|--------|-----|-------|------|----------|------|------|
| Chromic chloride | 10025-73-7 | 233-038-3 | -      | -   | X     | X    | KE-06017 | X    | X    |

| Component | CAS No | TSCA | TSCA Inventory notification - | DSL | NDSL | AICS | NZIoC | PICCS |
|-----------|--------|------|-------------------------------|-----|------|------|-------|-------|
|-----------|--------|------|-------------------------------|-----|------|------|-------|-------|

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

|                  |            |   |                        |   |   |   |   |   |
|------------------|------------|---|------------------------|---|---|---|---|---|
|                  |            |   | <b>Active-Inactive</b> |   |   |   |   |   |
| Chromic chloride | 10025-73-7 | X | ACTIVE                 | X | - | X | X | X |

**Legend:** X - Listed '-' - Not Listed

**KECL** - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

**Authorisation/Restrictions according to EU REACH** Not applicable

| Component        | CAS No     | REACH (1907/2006) - Annex XIV - Substances Subject to Authorization | REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances | REACH Regulation (EC 1907/2006) article 59 - Candidate List of Substances of Very High Concern (SVHC) |
|------------------|------------|---|---|---|
| Chromic chloride | 10025-73-7 | -   | -   | -   |

**Seveso III Directive (2012/18/EC)**

| Component        | CAS No     | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements |
|------------------|------------|---|--|
| Chromic chloride | 10025-73-7 | Not applicable  | Not applicable   |

**Regulation (EC) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of dangerous chemicals**

Not applicable

**Contains component(s) that meet a 'definition' of per & poly fluoroalkyl substance (PFAS)?**

Not applicable

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work .

## National Regulations

**UK** - Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment

**WGK Classification** See table for values

| Component        | Germany - Water Classification (AwSV) | Germany - TA-Luft Class |
|------------------|---------------------------------------|-------------------------|
| Chromic chloride | WGK1                                  |                         |

## 15.2. Chemical safety assessment

A Chemical Safety Assessment/Report (CSA/CSR) has not been conducted

## SECTION 16: OTHER INFORMATION

### Full text of H-Statements referred to under sections 2 and 3

H290 - May be corrosive to metals

H302 - Harmful if swallowed

H317 - May cause an allergic skin reaction

H411 - Toxic to aquatic life with long lasting effects

### Legend

# SAFETY DATA SHEET

Chromium(III) chloride, anhydrous

Revision Date 26-Jan-2024

**CAS** - Chemical Abstracts Service

**EINECS/ELINCS** - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**IECSC** - Chinese Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**ENCS** - Japanese Existing and New Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**NZIoC** - New Zealand Inventory of Chemicals

**WEL** - Workplace Exposure Limit

**ACGIH** - American Conference of Governmental Industrial Hygienists

**DNEL** - Derived No Effect Level

**RPE** - Respiratory Protective Equipment

**LC50** - Lethal Concentration 50%

**NOEC** - No Observed Effect Concentration

**PBT** - Persistent, Bioaccumulative, Toxic

**TWA** - Time Weighted Average

**IARC** - International Agency for Research on Cancer

Predicted No Effect Concentration (PNEC)

**LD50** - Lethal Dose 50%

**EC50** - Effective Concentration 50%

**POW** - Partition coefficient Octanol:Water

**vPvB** - very Persistent, very Bioaccumulative

**ADR** - European Agreement Concerning the International Carriage of Dangerous Goods by Road

**IMO/MDG** - International Maritime Organization/International Maritime Dangerous Goods Code

**OECD** - Organisation for Economic Co-operation and Development

**BCF** - Bioconcentration factor

**Key literature references and sources for data**

<https://echa.europa.eu/information-on-chemicals>

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

**ICAO/IATA** - International Civil Aviation Organization/International Air Transport Association

**MARPOL** - International Convention for the Prevention of Pollution from Ships

**ATE** - Acute Toxicity Estimate

**VOC** - (Volatile Organic Compound)

## Training Advice

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

Use of personal protective equipment, covering appropriate selection, compatibility, breakthrough thresholds, care, maintenance, fit and standards.

First aid for chemical exposure, including the use of eye wash and safety showers.

**Prepared By**

Health, Safety and Environmental Department

**Creation Date**

14-May-2010

**Revision Date**

26-Jan-2024

**Revision Summary**

New emergency telephone response service provider.

**This safety data sheet complies with Regulation UK SI 2019/758 and UK SI 2020/1577 as amended.**

## Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet**