

Rotor, adapters & Optional Accessories for CT-6E




Rotors

Rotor (Part No.)	Adapter pcs/box	Max. (rpm)	Max. RCF (Xg)	Max no of samples	Applicable tube	
T5SS (90585900)	5*7S	4	4800	4170	28	4/5mL blood tube
	10*7S	4	4800	4170	28	7mL blood tube
	15*7S	4	4800	4170	28	10mL blood tube
	10*7S	4	3000	1630	28	10mL glass tube
	15*7S	4	3000	1630	28	15mL glass tube
	50*1S	4	4800	3990	4	50mL glass tube
	15TC*4S	4	4800	4170	16	15mL culture tube (V bottom)
	15TC*3S	4	4800	4170	12	15mL culture tube (V bottom)
	50TC*1S	4	4800	4170	4	50mL culture tube (V bottom)
	160*1S	4	4800	4170	4	160mL bottle
	—	—	4800	4170	—	250mL bottle (Nalgene)
T4SS (90584100)	5*8S	4	4000	2380	32	4/5mL blood tube
	10*8S	4	4000	2860	32	7mL blood tube
	15*8S	4	4000	2860	32	10mL blood tube
	15*10S	4	4000	2860	40	10mL blood tube
	10*8S	4	3000	2320	32	10mL glass tube
	15*8S	4	3000	2320	32	15mL glass tube
	15*10S	4	3000	1610	40	15mL glass tube
	15TC*6S	4	4000	2900	24	15mL culture tube (V bottom)
	50TC*2S	4	4000	2860	8	50mL culture tube (V bottom)
T5SS3 (90591000)	2*48S	2	4800	3120	96	1.5/2.0mL micro tube
	—	—	4800	3120	2	96 well microplate
	—	—	4800	3120	2	96 deep well plate
	2*96MP	2	4800	3120	2	PCR 96 well plate
RT2S2 (90583000)	—	—	2000	700	4	96 well microplate
	15TA50T	2	6000	4630	8	15mL culture tube (V bottom)
T6AP (90591100)	—	—	6000	4770	8	50mL culture tube (V bottom)

Optional Accessories

Product	Model	Product No.	Remarks
Shield cap	S410252A	222480	4 pcs/set (T5SS/round bucket)
Round bucket	S410523A	222490	4 pcs/set (T4SS)
10mL glass tube	464977A	222500	50 pcs
15mL glass tube	464976A	222510	50 pcs
50mL glass tube	473866A	222520	10 pcs
160mL PP bottle	S306847A	222530	10 pcs
160mL PC bottle	S306848A	222540	10 pcs
160mL bottle cap	S410008A	222550	2 pcs



TOKYO HEADQUARTERS

Nihonbashi Hon-cho Bldg.,
3-3-4, Hon-cho Nihonbashi, Chuo-ku,
Tokyo 103-0023 Japan
TEL: 03-5201-6421
FAX: 03-3245-1225
e-mail: info.eyela@eyela.co.jp

US OFFICE

1190 Mt. View-Aviso Rd. Suite T
Sunnyvale, California 94089
TEL: 408-747-1543
FAX: 408-747-1523
e-mail: info@eyelusa.com

Tokyo Rikakikai Co., Ltd.

Singapore Liaison office
51 Anson Rd. #07-59
Anson Center, Singapore 079904
TEL: 6225-1306
FAX: 6225-1841
e-mail: eyelasingapore@singnet.com.sg

EYELA CO., LTD

Shanghai Sales
Rm 2209 Summit Center,
No.1088 West Yanan Rd.
TEL: 021-6207-6878
FAX: 021-6207-6920
e-mail: j-he@eyela.co.jp

Beijin Office
TEL: 010-8426-4221
FAX: 010-8426-4233
e-mail: fang-xiaodong@sohu.com

Guangzhou Branch
TEL: 020-3887-1910
FAX: 020-3887-1930
e-mail: j-tian@eyela.co.jp

