


AERATION UNIT Model MAU-1/2

*Powerful and compact compressor comes with the cabinet.
Specific air flow valve uniformly control air output rate.*



Combination :

- MAU-1C with upto 500ml capacity fermentor
- MAU-1D with upto 2.5L capacity fermentor
- MAU-2C with upto 5L capacity fermentor
- MAU-2D with upto 10L capacity fermentor

Specifications

Model	MAU-1C	MAU-1D	MAU-2C	MAU-2D
Air flow control range	: 50~500Nml/min.	200~2,500Nml/min.	0.5~5NL/min.	1~10NL/min.
Air compressor	: Max. 3.5NL/min.		Max. 10NL/min.	
Max. pressure	: 1.0kgf/cm ²		1.5kgf/cm ²	
Air filter	: Millipore filter FG-50			
Delivery tube	: Silicone tube (option)			
Cabinet dimensions	: 136Wx306Dx236H (mm)		166Wx356Dx301H (mm)	
Weight	: 4.5kg		7.5kg	
Power consumption	: 25VA		60VA	
Electricity	: 100VAC, 115VAC or 220/240VAC, 50/60Hz, 1 phase			
Replacement Parts	: Air filter, Delivery tube			

TOKYO RIKAKIKAI CO., LTD.

Tokyo, Japan

Tel. 81/3-3245-0481

Fax. 81/3-3241-0177