

EYELA

MAGNETIC STIRRER WITH BATH

MODEL;PS-100



PS-100

*Specifications

Speed Control:	100~1000rpm
Temp. Control Scale:	0~100°C
Temp. Control Range:	Room Temp.+5°C~95°C
Accuracy:	+-.0.1°C
Motor:	Induction type, 25W
Temp. Controller:	Thermistor Propotional Difference
Water Bath:	O250x110H 5L
Heating Element:	1kW
Power Source:	AC100V 50/60Hz 1.1kVA
Overall dimensions:	O250x220H
Standard Attachment:	Stirring Bar 1 pce Shelf 1 pce

TOKYO RIKAKIKAI CO.,LTD.