

Temperature Baths for Viscometers

Viscomate VM-150 III



Feature

- Small installation footprint, high accuracy circulating temperature bath
- Liquid contact sections of highly durable stainless steel (SUS304)
- Excessive load, overheat protection standard equipment
- Digital temperature setting and display

Specifications

Working temperature range	-10~+100°C
Temperature regulator accuracy	±0.05°C
External dimensions	W200 X D390 X H570mm
Weight	approx. 26 kg
Circulation pump	approx. 4L/min, lift head 1.4m
Maximum power consumption	approx. 1.0kVA
Power	AC100V ±10% 50/60Hz