



	NTT-20S	NTT-20G
Temp. control range	RT+5°C~ 95°C	RT+5°C~ 180°C
Temp. control accuracy	±0.05°C~	±0.1°C~
Temp. controller	PID contact less zero crossover	
Temp. setting/display	Membran key, increment 0.1°C(20G 1°C), Digital display	
Safety functions	self-diagnosis(upper/lower temp. abnormality detection, heater fault empty heating prevention, independent overheat protection, circuit protector	
Aux. functions	Timer function, alarm buzzer, auto tuning	
		Terminals for external sensor, recorder, alarm output
Heater	1kW (for water and oil)	
Bath agitation	Water jet agitation	
	Max. delivery 5L/min, Max. lift 1.2L	
Sensor	Themister	Pt100Ω
Ambient temp. for use	5~35°C	
Overall dimensions	194W x 306D x 359H	
Weight	Approx. 8kgs	
Power consumption	11A, 1.1kVA	
Rater voltage	1phase AC100V, 50/60Hz.	